



Sony Biotechnology Inc.

SA3800 Spectral Analyzer

The Sony SA3800 spectral cell analyzer incorporates advanced electronics and patented optical technologies to provide excellent analytical data. Its many automated features deliver true workflow simplicity.

Spectral technology in the SA3800 optimizes sensitivity and enhances dim signal detection by collecting photons at wavelengths from 420 nm to 800 nm. It also simplifies multicolor application design, workflow, and analysis for experienced and novice users.

Automation is present across all operations, from instrument startup, to quality control, acquisition, analysis, and system maintenance. Software wizards guide users through procedures across the workflow, making the system approachable for flow cytometry users of all levels.

Novel global standardization automatically sets the system to a master specification with a single click. This capability minimizes instrument variability from day to day and across multiple instruments for greater reliability.

The system features easy-to-use software and automation that speeds tasks across the workflow from setup to shutdown. Automated alignment, Experimental Designer, and Spectral Library are just a few of the the features that simplify operation. Each system includes fully integrated FCS Express™ software to offer the highest flexibility in analysis.

Spectral technology, easy-to-use software, and automation across the workflow make the SA3800 an excellent fit for facilities looking for a capable and easy-to-use analyzer for cellular acquisition and analysis.



- Uses spectral analysis technology to optimize sensitivity while simplifying application design and workflow.
- Unique 3D plate motion technology enables the plate (and sample) to move in horizontal and vertical directions to eliminate the need for cumbersome fluidics and their associated sample-to-sample cross-contamination. In addition, built-in, adjustable mixing mechanisms help maintain sample heterogeneity and optimal CVs between samples.
- The SA3800 features automation across the workflow, including automated alignment and a software calibration wizard to simplify operation and improve the reliability of results.



SA3800 System Overview


The SA3800 spectral analyzer improves sensitivity and simplifies application design, workflow, and analysis over conventional flow cytometers. This is achieved using spectral analysis technology, automation throughout the system, advanced electronics, and patented optical technologies, unique to Sony systems.

Supports up to 4 lasers including the 488-nm (blue) excitation laser that is standard in all systems. The 405-nm (violet), 561-nm (yellow-green), and 638-nm (red) lasers are available as options. 63.5 cm **The Flowpoint™ detection system** precisely tracks the core stream shape and position in the flow cell as well as the cross-sectional position of each passing particle to provide highly reliable measurements. This patented technology visualizes core stability and enables the highest resolution. Emitted light is collected by a 32-channel PMT producing a spectral fingerprint from the collected 32 data points of detected signal. A microlens array maximizes photon capture by A unique prism collection system delivers light refocusing light from the prism onto the 32through 10 consecutive prisms, to produce an channel PMT (photomultiplier tube) array optimal signal with minimal loss of light. 32-channel PMT Microlens array

Spectral Analysis Technology

Spectral analysis technology is the foundation of the SA3800 system. Spectral flow cytometry streamlines workflow and yields better data by collecting the entire visual spectrum of light. In conventional systems, overlapping fluorescence is subtracted using color compensation, so less light is collected. Instead, spectral flow cytometers sum the fluorescence together and then use unmixing to mathematically separate the colors. This powerful capability also simplifies workflow, including panel design, and improves visualization of autofluorescence.


A unique prism collection system delivers emitted light to a 32-channel PMT. This produces 34 data points of signal detection for fluorescence and bright autofluorescence to achieve accurate visualizations of fluorescent populations. Researchers can see the complete spectral fingerprint of each fluorochrome from 420 nm to 800 nm.



Visualization of fluorescent cell populations using the analyzer's unique prism collection system.

Spectral Unmixing

A powerful capability of spectral technology is unmixing. This technique allows researchers to separate fluorophores into pure signals that measure the quantity of each fluorophore at each pixel to more accurately measure data for analysis.




Spectral unmixing separates each spectral fingerprint for complete and optimal visualization of fluorochromes.

Spectral unmixing separates each spectral fingerprint to better visualize each fluorochrome marker. Unlike conventional filtering in which overlapping signals are lost, spectral unmixing captures the photons emitted from 420 nm to 800 nm. In doing so it enhances dim signal detection for better visualization of rare populations, fluorescent proteins, and fluorochromes excited by multiple lasers.

This also lets researchers separate the spectra of fluorophores masked by autofluorescence by extracting it from the signal and creating a unique autofluorescence parameter. With a clear signal for each color channel unaffected by overlapping signals (spillover) and autofluorescence, spectral analysis yields better, unbiased data for analysis.

Spectral analysis reduces false positives and delivers more accurate analysis over conventional flow cytometry. Mouse splenocytes were stained with CD4 APC and anti-IL-2 PE. **A**. In this conventional density plot, it is unclear if the light-blue region is a dim PE, weak double positive, or nonspecific binding.

B. Using spectral analysis, the spectral data of each region is compared against the contents of the Spectral Library to unmix the sample. This reveals the light-blue region is high autofluorescence. Representative data collected on the Sony SP6800.



System Automation

Automation is present across the workflow to simplify operation and ensure accurate results. The system supports a wide variety of standard and deep-well plates, and 12 x 75mm 5-mL tubes in the tube loader.

Novel 3D AutoSampler Technology

The novel 3D AutoSampler uses a fixed probe and moves the plate in horizontal and vertical directions to minimize sample-to-sample cross-contamination and speed up cleaning. Sensors on the probe enable the system to accurately move the sample to the probe, calculate height of the tube or plate, and automatically recoil the probe if it touches the container or base plate surface. This unique design is aimed at reducing cross-contamination and clogs, and allows for automatic recovery when flow problems are encountered. It also enables faster cleaning between samples.



Plate motion uses a fixed probe and moves the plate (and sample) in horizontal and vertical directions.







The SA3800 supports standard height and 96- and 384-well plates with round, flat, v-, and conical shapes in addition to 96-well half deep and bottom plates.








The tube loader supports 5-mL 12 x 75-mm tubes.


Mixing Automation and Cooling Plate

The 3D AutoSampler mixing function ensures consistent sampling throughout the acquisition of 96- or 384-well sample plates. Software optimizes settings for each plate type. Importantly, the mixing function maintains the integrity and heterogeneity of samples, ensuring that all particles are properly suspended for consistent results.

A cooling plate at the surface of the 3D base further reduces variability and inhibits sample degradation over time.



This figure illustrates how cell heterogeneity is maintained using the system's 3D mixing function. Without mixing, larger cells settle to form heterogeneous samples, as shown in the lower illustration. Graph data describes how the sample integrity is maintained throughout plate acquisition of 3-µm and 10-µm beads from wells 1–96 to deliver more complete, consistent results.



The graph describes the impact of the cooling plate on sample temperature during an 80-minute plate acquisition.

Software

The SA3800 software is easy to learn and use. It guides researchers from setup to panel design, acquisition, analysis, and shutdown.

User preferences allow operators and administrators to specify options for general instrument operation and how experiments are set up. These preferences promote ease of use and facilitate unattended operation procedures.

System Startup

At startup, Align Check and Performance QC wizards check instrument calibrations, using beads to ensure the instrument is operating optimally. On-screen instructions guide the user through procedures, then display progress and report results. The performance report displays MESF, Q, and B values to describe real-time fluorescence detection performance. If desired, Align Check and Performance QC reports can be displayed in historical context.

Standardization

Standardization mode sets the system to a master specification to minimize variability in a single instrument or among instruments located across sites. This unique capability allows experiments performed on any Sony cell analyzer to produce highly reliable, accurate, and reproducible results. This function also eliminates setup subjectivity for collaborative or long-term studies conducted by different operators or experience levels across sites.



When engaged, standardization mode sets the SA3800 to predetermined global settings that minimize instrument variability.

Spectral Library

The Spectral Library lets users create a personal reagent library that simplifies experiment creation and saves time. An Acquisition Wizard assists users with stepby-step instructions to acquire and analyze single-positive controls for the Spectral Library. Once acquired, the spectral reference for that reagent, including the spectral index, is available for future use. Information from the Spectral Library is available to users with a simple click, improving accuracy and streamlining panel design.



Experiment Creation

Experiments can be created using a template, an existing experiment, a singlestained, or a multicolor assay, in the Create Experiment window. Users can point and click to select (or edit) an existing experiment and can easily select templates for wells or plates when creating a new experiment. A setup Assay Wizard guides users through the creation of a single-stained or multicolor assay, simplifying experiment creation.



Acquisition Functions and Analysis

All acquisition functions, including instrument settings, are controlled from the Acquisition window. Worksheet tools let users choose how the data is displayed (such as plot types), and customize for their analysis needs. Plots and statistics provide real-time information during acquisition. To increase sample flow to the cuvette, a variable booster lets users set acquisition speeds from low (33 µL/min) normal, or high (250 µL/min) offering flexibility.



Walk-Away Event Checker


To enable walk-away operation, the Event Checker monitors stream integrity to identify air bubbles or low event rates that might be produced by clogs. When an issue is identified, the Event Checker preferences can be set to either clear the clog and stop, or clear the clog and continue to the next sample. A tube/well map shows which tube or well had the issue.

Autofluorescence Finder


In conventional flow cytometry, cellular autofluorescence produced by pyridine (NAD/NADH), flavin (FMN, FAD), and other intracellular oxidative reactions can cause fluorescence signal contamination in, for example, the FITC, PE, and Pacific Blue™ channels. Using the Autofluorescence Finder with spectral technology, autofluorescent spectral fingerprints can be easily found and subtracted to allow researchers to see the true fluorescent population.

Spectral Overlays

Comparative spectral curves can be created by using the Spectral Overlays function, to help identify unique populations between individual tubes/wells in the same or different experiments.







System Shutdown

On shutdown, the SA3800 software guides the user though pre-shutdown cleaning and shuts down the instrument automatically. Software wizards are also available to guide users through Bleach Cleaning and Rinse procedures.

FCS Express from De Novo Software

Spectral analyzers from Sony include FCS Express, from De Novo Software. FCS Express offers a range of new analysis tools from live gating to batch analysis. Native support for Sony spectral data files enables sophisticated data transformations and visualizations such as spectral overlays, tSNE, Spade, and heat maps.

Live Gating





Plate and High Content Analysis

Import plate or high content screening data and transform it into heat maps, overlays with radius visualization, and parameter statistics.




Flexible Analysis

Flexible data analysis can include integrated spreadsheets, custom calculations, charting, and regression analysis.

	1	c	5	1
18	FITC	PE	PerCP-Cy5.5	PE-Cy7
: M1	56.58	44.16	15.18	22.08
: M2	1818.84	1718.1	506.46	311.88
- M3	14735.64	14209.86	4364.94	2644.08
• M4	41860.92	40173.18	14921.94	• 10414.86
• M5	105084.24	101125.02	54574.86	46011.96
* M6	202089.95	204531.17	195759.91	• 198312.21


Expanded Data Visualization

Spectral, FCS, and CSV files are read directly by the software to enable more in-depth visualization of data.



Batch and Report

Batch analysis lets users process any number of samples quickly and easily with one click. To support presentation, data can be exported directly to PowerPoint, PDF, and Excel.




Sample Data

Using spectral unmixing, the SA3800 allows maximum flexibility for fluorochrome selection. Results are comparable to those obtained using conventional systems.

10-color Phenotyping of Human TBNK and Mature Erythrocytes using 3 Lasers

Alexa Fluor[®] 700


HLA-DR




2138125

Activated T cells

Figure A shows single positive control overlays of the spectral curves for CD45, CD3, CD4, and CD8 from a human TBNK sample. When compared to the TBNK sample (Figure B), the single positive control overlays match. CD4 Alexa Fluor® 647 and CD8 Alexa Fluor® 532 are used to distinguish the presence of T-helper cells and T-cytotoxic cells in the sample.





10

10'


10³ -10³

Intensity_H

Sample Data


Comparative Data on Human Lymphocytes

A four-color comparative analysis of fixed human lymphocytes was done to evaluate resolution and sensitivity. The same sample was used on the same day on four different instruments. Each instrument had the recommended cleaning and instrument QC performed before running the sample. When compared with competitors, the Sony SA3800 and SP6800 spectral analyzer results showed better separation between the negative and positive populations. In addition, a comparison was performed of the removal of autofluorescence. Each instrument showed autofluroescent populations stemming from the negative population. The SA3800 was able to identify the population, create it as its own color using the Autofluorescence Finder, and subtract the populations to allow the user to see the true populations.



4-Color Comparative Data on Human Lymphocytes

Comparative Data with Autofluorescence Removal



Comparative Data on Human Lymphocytes


Marker	Fluorochrome	Purpose	Cat. No.
CD45	PerCP-Cy™5.5	Leukocytes	2120135
CD16	PE	Activated cells	2110035
CD3	Alexa Fluor [®] 488	T cells	2102075
CD8	APC-Cy7	Effector T cells	2105075




Sample Data

Using 10 Colors with 405-nm and 488-nm Lasers to find Multiple Populations

In an experiment using fixed human blood, B cells (CD19), macrophages (CD33), basophils (CD123), NK cells (NKp46), and helper T cells (CD45/CD3/CD4) were identified. Also, a smaller subset population of NKT cells was analyzed using CD57. CD57 shows expression mostly on negative CD4 cells, since it is associated with CD8 killer T cells.






10 Colors to find Multiple Populations

Marker	Fluorochrome	Purpose	Cat. No.
CD45	Alexa Fluor® 532	Leukocytes	N/A
CD3	PerCP-Cy5.5	T cells	2186675
CD4	PE Alexa Fluor® 700	Helper T cells	2187125
CD19	PE/Dazzle™ 594	B cells	2111255
CD33	PE-Cy7	Macrophages	N/A
CD123	BV421	Basophils	2130085
CD11c	PE	Dendritic cells	2108025
CD335 (NKp46)	BV650	NK cells	2259635
CD57	FITC	NKT-cell subset	2211530
HLA-DR	Pacific Blue™	Macrophages, basophils	2138115

Spectral Fingerprinting Using a TBNK Panel

Cytometrists use many forms of controls to determine cell populations from a heterogeneous mixture of cells. In this example, a 6-color TBNK panel was analyzed completely by spectral fingerprints. By comparing the single positive spectral fingerprints, the need for an FMO control to identify a population is decreased. In addition, since the autofluorescence from the unstained cells is treated as a separate color, it is subtracted to be able to identify the cellular lineage easily.

Single Positive Control Spectral Graphs



Specifications

	ltem	Specifications
	Excitation options	405 nm, 488 nm, 561 nm, 638 nm
	Detection signals	FSC, SSC, FL 32chPMT (500 to 800 nm), Violet x2ch PMT (420-500 nm)*
	Pulse parameter	Area, Height, Width (all channels)
Optics/performance	Signal resolution	Height 20 bit, Area 32 bit/Sampling frequency: 50 MHz
	Fluorescence sensitivity	FITC: 120 MESF; PE: 70 MESF (nominal)
	Fluorescence resolution	CV <3% for the singlet peak of propidium iodide-stained CEN
	Detectable cell size	0.5 μm to 40 μm
	Event rate	20,000 eps (max)
	Sample volume rate	33 μL/min-240 μL/min
Fluidics	Cleaning	Auto probe cleaning, priming, shutdown cleaning, cuvette backflushing
	Waste/DI water tank	Both tanks are 2 L, 10 L available as an option
Single loader	Single tube	Falcon 5 mL (12 x 75-mm) polystryene/polypropylene
	Sample well	96-well plate: standard height flat/V/U, 96 half deep, 96 deep, 384 standard flat
	Sample tube	Tube rack: 24 Falcon 5-mL (12 x 75-mm) polystyrene/polypropylene
3D AutoSampler	Sample volume	Tube: minimum 100 μL/maximum 2,000 μL 96-well-plate: standard height: 55 μL–200 μL 384-well-plate: standard height: 40 μL–75 μL
	Carryover	<0.1% (under high speed/normal mode)
	Measurement speed	96-well plate in 25 minutes (Acquisition time per well: 2 seconds)
	Reagent stability	Mixing function, sample cooling block (passive)
	Dimensions	W: 66 cm x H: 67.4 cm x D: 63.5 cm (SA3800 main body)
Ancillary	Weight	95 kg (AutoSampler model, does not include external tank holder)
	Power consumption	350 W Max

3D AutoSampler model

Model	No. of lasers	Laser wavelengths (nm)
LE-SA3800AA	1	488
LE-SA3800BA	2	488, 638
LE-SA3800CA	2	488, 405
LE-SA3800DA	2	488, 561
LE-SA3800EA	3	488, 405, 638
LE-SA3800GA	3	488, 405, 561
LE-SA3800FA	4	488, 405, 638, 561

Single loader model

Model	No. of lasers	Laser wavelengths (nm)
LE-SA3800AS	1	488
LE-SA3800BS	2	488, 638
LE-SA3800CS	2	488, 405
LE-SA3800DS	2	488, 561
LE-SA3800ES	3	488, 405, 638
LE-SA3800GS	3	488, 405, 561
LE-SA3800FS	4	488, 405, 638, 561



North America/International

1730 North First Street San Jose, CA 95112 U.S.A. Voice: +1 800-275-5963 FAX: +1 408-352-4130 sales@sonybiotechnology.com http://www.sonybiotechnology.com

Japan

4-14-1, Asahi-cho, Atsugi-shi, Kanagawa, 243-0014 Japan Tel: +81 120-677-010 Fax: +81 120-388-060 sales_Japan@sonybiotechnology.com http://www.sony.co.jp/LS



The Heights, Brooklands. Weybridge, Surrey, KT13 0XW, UK sales_EU@sonybiotechnology.com

©2018 Sony Biotechnology Inc. All rights reserved. Sony, the Sony logo, and FlowPoint are trademarks of Sony Corporation. Pacific Blue™ is a trademark of, and Alexa Fluor® is a registered trademark of and licensed under patents assigned to Life Technologies Corporation. PE/Dazzle is a trademark of BioLegend. Cy is a trademark of GE Healthcare. De Novo and FCS Express are trademarks of De Novo Software. All other trademarks are property of their respective owners. For non-clinical research use only. Not for use in diagnostic or therapeutic procedures, or for any other clinical purpose. The SA3800 Spectral Analyzer is classified as a Class 1 Laser Product. Specifications subject to change without notice.